



Product Specification

Product : Acetone
Catalog# : GC9802
CAS : 67-64-1
Formula : C₃H₆O
MW : 58.08
Grade : GC

TEST	SPECIFICATION
Description	Colorless liquid
Purity (By GC)	≥ 99.8%
Identification by IR	To pass the test
Color (APHA)	≤ 10
Residue after evaporation	≤ 10ppm
GC-ECD	10 ng/l
Water	≤ 0.5%

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

