



Product Specification

Product : Benzene
Catalog# : GC9804
CAS : 71-43-2
Formula : C₆H₆
MW : 78.11
Grade : GC

TEST	SPECIFICATION
Description	Colorless clear liquid.
Purity (By GC)	≥ 99.7%
Identification by IR	To pass the test
Color (APHA)	≤ 10
Residue after evaporation	≤ 10ppm
Substance darkened by sulfuric acid	To pass the test
Water	≤ 0.05%
Sulfur compounds (as S)	≤ 0.005%
Thiophene	To pass the test

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

