



Product Specification

Product : Dichloromethane
Catalog# : GC9808
CAS : 75-09-2
Formula : CH₂Cl₂
MW : 84.93
Grade : GC

TEST	SPECIFICATION
Description	Clear colorless liquid free of suspended matter.
Assay (GC-FID)	≥ 99.8%
GC-ECD	≤ 10 ppt
Color (APHA)	≤ 10
Residue after evaporation	≤ 10 ppm
Stabilizer (As Amylenes)	Between 100-200 ppm
Chloride (Cl)	≤ 0.0010%
Water	≤ 0.02%

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

