



Product Specification

Product : Dimethyl sulfoxide
Catalog# : GC9810
CAS : 67-68-5
Formula : (CH₃)₂SO
MW : 78.13
Grade : GC

TEST	SPECIFICATION
Description	Colorless liquid.
Purity (By GC)	≥ 99.8%
Identification by IR	To pass the test
Residue after evaporation	≤ 0.002%
GC-ECD	≤ 10 ng/l
Water	≤ 0.2%
UV absorbance 350 nm	≤ 0.01
UV absorbance 300 nm	≤ 0.10
UV absorbance 280 nm	≤ 0.30
UV absorbance 270 nm	≤ 0.70

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

