



Product Specification

Product : Ethyl acetate
Catalog# : GC9811
CAS : 141-78-6
Formula : C₄H₈O₂
MW : 88.11
Grade : GC

TEST	SPECIFICATION
Description	Colorless liquid
Purity (By GC)	≥ 99.8%
Identification by IR	To pass the test
Color (APHA)	≤ 10
Residue after evaporation	≤ 5 ppm
GC-ECD	≤ 10 ppt
Substance darkened by sulfuric acid	≤ 0.0009 meq/g
Water	≤ 0.02%

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

