



Product Specification

Product : n-Hexane
Catalog# : GC9812
CAS : 110-54-3
Formula : C₆H₁₄
MW : 86.18
Grade : GC

TEST	SPECIFICATION
Description	Colorless liquid
Infrared spectrum	To conform the structure
Purity (By GC)	
As C6 hydrocarbons	≥ 99.5%
As n-Hexane	≥ 95.0%
Residue after evaporation	≤ 0.0001%
GC-ECD	≤ 10 ppt
GC-FID	≤ 5 ppb
Water	≤ 0.01%
Thiophene	To pass the test
Sulfur compounds (as Sulfur)	≤ 0.005%

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

