



Product Specification

Product : Isopropyl alcohol
Catalog# : GC9814
CAS : 67-63-0
Formula : C₃H₈O
MW : 60.10
Grade : GC

TEST	SPECIFICATION
Description	Colorless liquid
Purity (By GC)	≥ 99.8%
Identification by IR	To pass the test
Color (APHA)	≤ 10
Residue after evaporation	≤ 10ppm
GC-ECD	≤ 5ppt
Water	≤ 0.2%

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

