



Product Specification

Product : Methanol
Catalog# : GC9815
CAS : 67-56-1
Formula : CH₄O
MW : 32.04
Grade : GC

TEST	SPECIFICATION
Description	Clear Colorless liquid free of suspended matter
Purity (By GC)	≥ 99.8%
Identification by IR	To pass the test
Color (APHA)	≤ 10
Residue after evaporation	≤ 5 ppm
GC-ECD	≤ 5 ppt
Water	≤ 0.1%

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

