



# Product Specification

Product : N, N dimethylformamide  
Catalog# : GC9819  
CAS : 68-12-2  
Formula : C<sub>3</sub>H<sub>7</sub>NO  
MW : 73.09  
Grade : GC

TEST	SPECIFICATION
Appearance	Colorless liquid
Assay (by GC-FID)	≥ 99.80%
Color (APHA)	≤ 10
Residue after Evaporation	≤ 5ppm
GC-ECD	≤ 10 ng/l
Water	≤ 0.05 %
UV Absorbance @ 270 nm	1.00 AU
UV Absorbance @ 275 nm	0.30 AU
UV Absorbance @ 295 nm	0.10 AU
UV Absorbance @ 310 nm	0.05 AU
UV Absorbance @ 340 nm	0.01 AU
UV Absorbance @ 400 nm	0.01 AU

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

