



Product Specification

Product : Bromocresol Purple (Free Acid)
Catalog# : I10505
CAS : 115-40-2
Formula : $C_{21}H_{16}Br_2O_5S$
MW : 540.22

TEST	SPECIFICATION
Description	Grey to light purple or pink to peach to yellow-tan powder
Assay	≤ 94.0%
Melting point	241°C to 242°C
Clarity of solution	Passes test
Aqueous solubility	Practically insoluble
Visual Transition Interval	From pH 5.2 (Yellow) To pH 6.8 (Purple)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Tri Stains

