



Product Specification

Product : Methyl orange
Catalog# : I10538
CAS : 547-58-0
Formula : -
MW : -

TEST	SPECIFICATION
Description	Light to dark orange powder
Clarity of solution	To pass ACS test
Aqueous solubility	Soluble in 500 parts water; more soluble in hot water
Loss on drying at 110°C / 1 Hr	< 3.0%
Visual transition interval	From pH 3.2 (pink to light red) to pH 4.4 (yellow)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Tri Stains

