



Product Specification

Product : Methyl Red HCL
Catalog# : I10540
CAS : 63451-28-5
Formula : $C_{15}H_{15}N_3O_2 \cdot HCl$
MW : 305.77

TEST	SPECIFICATION
Appearance	Dark purple crystalline powder
Loss on drying	10% maximum
Clarity of alcohol Solution	Passes test
Visual Transition Interval	pH 4.2 (pink) to pH 6.2 (yellow)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

