



Product Specification

Product : Phenol Red Sodium Salt
Catalog# : I10556
CAS : 34487-61-1
Formula : C₁₉H₁₃NaO₅S
MW : 376.36

TEST	SPECIFICATION
Description	Dark-red to brownish-red powder
Identification (IR)	To pass test
Assay: (UV vis)	>90%
Loss on drying	< 5%
Clarity of solution	To pass Test
Lambda Max	423 nm in Methanol
Melting Point	285 ° C
Visual transition Interval	From pH 6.8 (yellow) to pH 8.2 (red)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Tri Stains

