



Product Specification

Product : Thymol blue sodium salt
Catalog# : I10570
CAS : 62625-21-2
Formula : $C_{27}H_{29}NaO_5S$
MW : 488.57

TEST	SPECIFICATION
Appearance	Red-Brown to brown of black powder
IR spectrum	Confirms the test
pH test	pH 1.2 (Red) to pH 28 (Yellow); pH 8.0 (Yellow) to pH 9.2 (Blue)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

