



Product Specification

Product Information:

Products:	Bromophenol Blue 0.04% solution
Catalog:	IS10760
Grade:	Indicator solution

Reagent Used:

Reagents	CAS #	Grade
Bromophenol blue sodium salt	62625-28-9	ACS
Water	7732-18-5	USP/ASTM

Specification:

Test	Specification
Appearance	Dark reddish-purple liquid
Suitability for Use	Yellow (pH < 3.0) – Green (pH 3.4 to 3.7) – Blue (pH >4.6)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Tri Stains

