



# Product Specification

## Product Information:

<b>Products:</b>	Bromocresol Green 0.1% (Aqueous Solution)
<b>Catalog:</b>	IS10766
<b>Grade:</b>	Indicator solution

## Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	ASTM/USP/ACS
Bromocresol Green, Sodium Salt	62625-32-5	ACS

## Specification:

Test	Specification
Appearance	Dark blue liquid
Suitability for Use	Blue (pH >5.4) - Green (pH 4.5)- Yellow (pH < 4.0)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Tri Stains

