



# Product Specification

## Product Information:

<b>Products:</b>	Methyl Red 0.1%
<b>Catalog:</b>	IS10790
<b>Grade:</b>	Indicator solution

## Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	ACS/ASTM/USP
Methyl Red, Sodium Salt	845-10-3	ACS

## Specification:

Test	Specification
Appearance	Orange to dark red liquid
Suitability for Use	Pink (pH 4.5) - Yellow (pH 6.2)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Tri Stains

