



Tri Stains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product Information:

Products:	Phenolphthalein 1% w/v in Ethanol
Catalog:	IS10806
Grade:	Indicator solution

Reagent Used:

Reagents	CAS #	Grade
Ethyl alcohol	64-17-5	Reagent
Methyl alcohol	67-56-1	Reagent
Phenolphthalein	77-09-8	ACS

Specification:

Test	Specification
Appearance	Colorless to slightly yellow liquid
Suitability for use	Colorless (pH 8.0) – Pink (PH 10.0)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Tri Stains

