



Product Specification

Product Information:

Products:	Thymolphthalein 0.05% in 90% Alcohol
Catalog:	IS10813
Grade:	Indicator solution

Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	ASTM/USP
Thymolphthalein	125-20-2	ACS
Isopropyl Alcohol	67-63-0	Reagent

Specification:

Test	Specification
Appearance	Colorless liquid
Suitability for use	Colorless (pH < 8.8) – Pale, greyish-blue (pH 9.4) - Blue (pH 9.9) - Intense Blue (pH > 10.5)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Tri Stains

