



Product Specification

Product Information:

Products:	Phenolphthalein 1.0% in 50% Isopropyl alcohol
Catalog:	IS10821
Grade:	Indicator solution

Reagent Used:

Reagents	CAS #	Grade
Isopropyl Alcohol	67-63-0	ACS/USP
Water	7732-18-5	USP/ASTM
Phenolphthalein	77-09-8	ACS

Specification:

Test	Specification
Appearance	Colorless to slightly pink liquid
Suitability for Use	Colorless (pH 8.0) - Pink (pH 10.0)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Tri Stains

