



Product Specification

Product Information:

Products:	Methyl Red 0.2%
Catalog:	IS10822
Grade:	Indicator solution

Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	ACS/ASTM/USP
Methyl Red, Sodium Salt	845-10-3	ACS

Specification:

Test	Specification
Appearance	orange to dark red liquid
Suitability for Use	Pink (pH 4.5) - Yellow (pH 6.2)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Tri Stains

