



Tri Stains

**DAWN SCIENTIFIC**

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

# Product Specification

## Product Information:

<b>Products:</b>	Crystal Violet Indicator, 0.2% in Glacial Acetic Acid
<b>Catalog:</b>	IS10828
<b>Grade:</b>	Indicator solution

## Reagent Used:

Reagents	CAS #	Grade
Acetic Acid	64-19-7	ACS
Gentian Violet	548-62-9	ACS

## Specification:

Test	Specification
Appearance	Dark blue to violet/purple liquid
Suitability for Use	Yellow (pH < 0.1) - Blue-Violet (pH > 2.0)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Tri Stains

