



Tri Stains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product Information:

Products:	Murexide Indicator 0.15%
Catalog:	IS10841
Grade:	Indicator solution

Reagent Used:

Reagents	CAS #	Grade
Ethylene glycol	107-21-1	ACS
Murexide	3051-09-0	ACS

Specification:

Test	Specification
Appearance	Dark red liquid
Suitability for use	Pink (Excess Ca Ions) – Purple (Excess EDTA Ions)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Tri Stains

