



Product Specification

Product Information:

| | |
|------------------|---|
| Products: | Phenolphthalein 1% in 60% Isopropyl alcohol |
| Catalog: | IS10845 |
| Grade: | Indicator solution |

Reagent Used:

| Reagents | CAS # | Grade |
|-------------------|-----------|----------|
| Isopropyl Alcohol | 67-63-0 | ACS/USP |
| Water | 7732-18-5 | USP/ASTM |
| Phenolphthalein | 77-09-8 | ACS |

Specification:

| Test | Specification |
|---------------------|-------------------------------------|
| Appearance | Colorless to slightly pink liquid |
| Suitability for Use | Colorless (pH 8.0) - Pink (pH 10.0) |

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Tri Stains

