



Product Specification

Product Information:

| | |
|------------------|---------------------------------------|
| Products: | Starch 1% (Preserved with Chloroform) |
| Catalog: | IS10848 |
| Grade: | Indicator solution |

Reagent Used:

| Reagents | CAS # | Grade |
|----------------|-----------|----------|
| Starch soluble | 9005-84-9 | ACS |
| Water | 7732-18-5 | USP/ASTM |
| Chloroform | 67-66-3 | ACS |

Specification:

| Test | Specification |
|---------------------|---|
| Appearance | Translucent liquid |
| Suitability for Use | Colorless (Iodine absent) – Blue (Iodine present) |

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

