

Product Specification

Product Information:

Products:	Ammonium chloride 10% solution
Catalog:	PS18009
Grade:	Percentage solution

Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	USP/ASTM
Ammonium chloride	12125-02-9	ACS

Specification:

Test	Specification
Appearance	Colorless liquid.
Assay	9.7 to 10.3%

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Verified By



Anjal Patel

