



Product Specification

Product Information:

Products:	Potassium carbonate 15% w/v
Catalog:	PS18083
Grade:	Percentage solution

Reagent Used:

Reagents	CAS #	Grade
Potassium carbonate	584-08-7	ACS
Water	7732-18-5	USP/ASTM

Specification:

Test	Specification
Appearance	Colorless liquid.
Assay	15.0 % \pm 0.5%

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Verified By

Anjal Patel

