

Product Specification

Product Information:

Products:	Potassium fluoride 50% solution
Catalog:	PS18094
Grade:	Percentage solution

Reagent Used:

Reagents	CAS #	Grade
Potassium fluoride	7789-23-3	Reagent
Water	7732-18-5	USP/ASTM

Specification:

Test	Specification
Appearance	Clear colorless liquid
Assay	50.0 % ± 1.0%

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Verified By



Anjal Patel

