



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Acid blue 29
Catalog# : S10002
CAS : 20460
CI number : 5850-35-1
Formula : $C_{22}H_{14}N_6O_9S_2Na_2$
MW : 616.50

TEST	SPECIFICATION
Appearance	Dark brown to navy blue or black powder
Identification	Conforms to structure and Standard
UV-VIS (0.015-0.024 g/l, H ₂ O)	
@599-605nm	NLT 15000
@388-394nm	NLT 3000
@308-314nm	NLT 8000
@230-236nm	NLT 9000
Solubility (0.1-0.14% in H ₂ O)	Opaque, dark blue solution

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

