



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Acid Violet 17
Catalog# : S10006
CAS : 4129-84-4
CI number : 42650
Formula : $C_{41}H_{44}N_3O_6S_2Na$
MW : 761.92

TEST	SPECIFICATION
Appearance	Dark Purple to dark green to green brown or bluish violet crystalline powder
Identification (IR)	Conforms to structure
Dye content	$\geq 48.0\%$
UV Vis: (c = 0.02 g/l)	
543 - 549nm	$\geq 34,000$
304 - 310nm	$\geq 9,000$
250 - 256nm	$\geq 9,000$
Solubility	
H2O	30 mg/ml
EGME	30 mg/ml
EtOH	6 mg/ml

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

