



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Alizarin Yellow R
Catalog# : S10018
CAS : 2243-76-7
CI number : 14030
Formula : $O_2NC_6H_4N=NC_6H_3-2-(OH)CO_2H$
MW : 287.23

TEST	SPECIFICATION
Appearance	Reddish Orange Powder
Solubility (0.1% IN Ethanol)	Clear Yellow Solution
Absorptivity In Ethanol	(A 1%/1cm) max-385nm-1004
Loss on drying (110 °C)	Max. 3%
Melting Point	>300 deg. C
Dye Content	Min. 75%

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

