



TriStains

**DAWN SCIENTIFIC**

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

## Product Specification

Product : Alizarin Yellow R Sodium Salt  
Catalog# : S10019  
CAS : 1718-34-9  
CI number : 14030  
Formula :  $C_{13}H_8N_3NaO_5$   
MW : 309.21

TEST	SPECIFICATION
Description	Reddish-orange powder
Solubility in water	Sparingly soluble
Melting Point	253-254°C
Dye content	96%
Solubility (0.1% in ethanol)	Clear yellow solution
Loss on drying (@ 110°C)	1.74%

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

