



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Aniline Blue W/S Certified
Catalog# : S10024
CAS : 28983-56-4
CI number : 42780
Formula : Not Applicable
MW : Not Applicable

TEST	SPECIFICATION
Appearance	Brownish-purple crystalline powder with metallic lustre
Solubility (in water)	Soluble
Lambda max.	594-610nm in water

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

