



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Auramine O
Catalog# : S10026
CAS : 2465-27-2
CI number : 41000
Formula : $C_{17}H_{22}N_3Cl$
MW : 303.84

TEST	SPECIFICATION
Appearance	Yellow to yellow green powder
Dye content	$\geq 80.0\%$
Absorption maxima @ 435nm	≥ 25000
Loss on drying	$\leq 5.0\%$

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

