



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Congo Red Certified
Catalog# : S10054
CAS : 573-58-0
CI number : 22120
Formula : $C_{32}H_{22}N_6Na_2O_6S_2$
MW : 696.67

TEST	SPECIFICATION
Appearance	Red to red-brown powder
Color Change	Blue-violet @ pH3.0 to red-orange @ pH 5.2
Loss on Drying	$\leq 5.0\%$
UV-Vis	490-500nm @ pH 7
Solution Test	0.1% in water: clear

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

