



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Eosin B Certified
Catalog# : S10066
CAS : 548-24-3
CI number : 45400
Formula : $C_{20}H_6N_2O_9Br_2$
MW : 624.08

TEST	SPECIFICATION
Appearance	Brownish-Green Powder
Lambda max	514 (395) nm (in Water plus 1% Na ₂ CO ₃)
Visual Transition Interval	pH 0.0 (colorless) to pH 2.4 (fluorescent pink)
Dye Content	85 % minimum

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

