



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Eosin B Spirit Soluble
Catalog# : S10067
CAS : 56360-46-4
CI number : 45400
Formula : $C_{20}H_8N_2O_9Br_2$
MW : 580.11

TEST	SPECIFICATION
Appearance	Yellow to Orange Powder
Identification (IR)	Conforms to standard
Absorptivity: 1% 1cm(in H ₂ O + 1ml, 2% Na ₂ CO ₃ at λ max. 515-524nm)	>1020
Solubility: 1mg/ml, MEOH (HEAT)	Clear Red solution
Visual Transition Interval	pH 1.4 (colorless) to pH 2.4 (fluorescent pink)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

