



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Eosin Y Certified
Catalog# : S10068
CAS : 17372-87-1
CI number : 45380
Formula : $C_{20}H_6Br_4Na_2O_5$
MW : 691.85

TEST	SPECIFICATION
Description	Dark red-brown (rust) powder
Identification (IR)	To conform to standard
Assay	> 88%
Solubility: (0.25%, in water)	Clear, yellow-orange to red solution
Loss on drying	< 5%
Certification	To be certified by the BIOLOGICAL STAINS COMMISSION

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

