



TriStains

**DAWN SCIENTIFIC**

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

## Product Specification

Product : Eriochrome Black T ACS  
Catalog# : S10070  
CAS : 1787-61-7  
CI number : 14645  
Formula :  $C_{20}H_{12}N_3NaO_7S$   
MW : 461.38

TEST	SPECIFICATION
Appearance	Black powder
Loss on drying	< 7%
Absorptivity 1%/1cm (At pH 10 Ammonia Buffer + 2ml 0.05 E.D.T.A. Solution) at $\lambda$ max. 623.5 nm	>320
Clarity of Solution	To pass test
Suitability as Complexometric Indicator	To pass test

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

