



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Fluorescein Diacetate
Catalog# : S10087
CAS : 596-09-8
CI number : 45330:1
Formula : $C_{24}H_{16}O_7$
MW : 416.39

TEST	SPECIFICATION
Description	Light yellow crystalline powder
Identification (IR)	Conforms to structure and standard
Assay (HPLC)	$\geq 97.0\%$
Melting Point	200°C – 205°C
Sulfated Ash	$\leq 0.1\%$

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

