



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Janus Green B Certified
Catalog# : S10102
CAS : 2869-83-2
CI number : 11050
Formula : $C_{30}H_{31}ClN_6$
MW : 511.07

TEST	SPECIFICATION
Description	Very dark purple or very dark green to almost black powder
Dye content	> 50%
Lambda max. (in 50% ethanol)	660 nm
Solubility	Soluble in water, slightly soluble in alcohol

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

