



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Leuco Crystal Violet
Catalog# : S10106
CAS : 603-48-5
CI number : Not applicable
Formula : $C_{25}H_{31}N_3$
MW : 373.53

TEST	SPECIFICATION
Appearance	White to light grey to very pale Lavender crystalline powder
Assay	98 % minimum
Lambda Max	261 nm in 0.1 N Hydrochloric acid in Methanol
Melting Point	175 - 180 ° C

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

