



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Methyl Blue
Catalog# : S10116
CAS : 28983-56-4
CI number : 42780
Formula : $C_{37}H_{27}N_3O_9S_3Na_2$
MW : 799.80

TEST	SPECIFICATION
Appearance	Brownish-Purple crystalline powder With metallic lustre
Solubility (in water)	Soluble
Lambda max	594-610 nm in water

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

