



TriStains

**DAWN SCIENTIFIC**

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

## Product Specification

Product : Methylene Blue Certified  
Catalog# : S10117  
CAS : 7220-79-3  
CI number : 52015  
Formula :  $C_{16}H_{18}ClN_3S \cdot 3H_2O$   
MW : 373.90

TEST	SPECIFICATION
Appearance	Dark greenish crystalline powder
Identification	To pass test
Assay (as is bases)	> 82%
Lambda max.	663-667nm in water, Ratio 1.21-1.7
Solubility	Soluble in water, ethanol, methyl Cellosolve and ethylene glycol.
Material to be certified by the Biological Stains Commission	

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

