



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Methylene Blue Anhydrous
Catalog# : S10118
CAS : 61-73-4
CI number : 52015
Formula : $C_{16}H_{18}ClN_3S$
MW : 319.85

TEST	SPECIFICATION
Appearance	Dark greenish crystalline powder
Identification	To pass test
Assay (Dried Bases)	98 - 103 %
Loss on drying (Dried at 100° C for 1Hr.)	3% maximum
Residue on ignition	1.2 % maximum
Arsenic	8 ppm maximum
Copper or Zinc	0.02 % maximum

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

