



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Methyl Violet 2B Certified
Catalog# : S10123
CAS : 8004-87-3
CI number : 42535
Formula : $C_{24}H_{26}N_3Cl$
MW : 393.96

TEST	SPECIFICATION
Appearance	Green Crystalline Powder
Assay	75% minimum
Lambda max	583-587 nm (in 50% Ethanol)
Visual Transition Interval	pH 0.1 (yellow) to pH 2.7 (pink)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

