



TriStains

**DAWN SCIENTIFIC**

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

## Product Specification

Product : Naphthol Green B ACS  
Catalog# : S10124  
CAS : 19381-50-1  
CI number : 10020  
Formula :  $C_{30}H_{15}FeN_3Na_3O_{15}S_3$   
MW : 878.450

TEST	SPECIFICATION
Description	Dark green powder
Absorptivity (A 1 cm 1% at 715 nm water, on dried substance)	105 – 115
Loss on drying (105°C)	≤ 10.0%
Solubility in water	Soluble in water

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

