



TriStains

**DAWN SCIENTIFIC**

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

## Product Specification

Product : Neutral Red Certified  
Catalog# : S10126  
CAS : 553-24-2  
CI number : 50040  
Formula :  $C_{15}H_{17}ClN_4$   
MW : 288.78

| TEST                                             | SPECIFICATION                           |
|--------------------------------------------------|-----------------------------------------|
| Appearance                                       | Very Dark-Green Powder                  |
| Melting Point                                    | ~290 ° C                                |
| Dye content                                      | > 50%                                   |
| Solubility 0.1% (50% Ethanol + 2 ml Acetic Acid) | Clear Red solution                      |
| Visual Transition Interval                       | From pH 6.8 (red)<br>To pH 8.0 (yellow) |

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

