



TriStains

**DAWN SCIENTIFIC**

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

## Product Specification

Product : Orange II Certified  
Catalog# : S10138  
CAS : 633-96-5  
CI number : 15510  
Formula :  $C_{16}H_{11}N_2NaO_4S$   
MW : 350.33

TEST	SPECIFICATION
Appearance	Orange-Brown Powder
Assay	87 % min
Lambda Max	484 - 487 nm

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

