



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Orange G Certified
Catalog# : S10140
CAS : 1936-15-8
CI number : 16230
Formula : $C_{16}H_{10}N_2Na_2O_7S_2$
MW : 452.38

TEST	SPECIFICATION
Appearance	Red-orange powder
Lambda max:	475 nm (in water)
Dye content	80% minimum
Solubility	Soluble in water and methyl Cellosolve Slightly soluble in ethanol

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

